



## esco HD-NaCl

### Sodium chloride in pharmacopoeia quality

- ✓ Ph. Eur., USP, JP
- ✓ Production in batches
- ✓ Analysis certificate

#### Field of application

esco HD-NaCl has been developed especially for use in hemodialysis. Other pharmaceutical applications are also possible.

#### Product properties

esco HD-NaCl is a particularly pure vacuum salt that is produced for demanding pharmacological applications. esco HD-NaCl is manufactured and tested in accordance with the requirements of the latest versions of the European, American and Japanese pharmacopoeias. Other pharmacopoeias are available upon request.

- ◆ High purity
- ◆ Of natural origin
- ◆ Low content of by-minerals



# esco HD-NaCl

## Sodium chloride in pharmacopoeia quality

### Product range

Packaging size: 25 kg bag  
(40 x 25 kg/  
IPPC pallet)  
50 kg bag  
(21 x 50 kg/  
IPPC pallet)  
1,000 kg  
single-use FIBC/  
Big Bag  
Bulk

Type of packaging: 25 kg, 50 kg:  
Pharma bag  
(multi-layer  
paper bag  
(with PE coating)  
or pure PE bag)  
1,000 kg  
single-use FIBC/  
Big Bag  
with inliner

GTIN: 25 kg:  
4003885449175  
50 kg:  
4003885449168

We would be happy to provide you  
with a sample.

For more information on this and  
other products, please do not hesitate  
to contact us.

### Analysis

esco HD-NaCl is produced in individual batches and then analyzed. With each delivery, an individual analysis certificate is provided. Accompanying HACCP analyses in accordance with the principles of the Codex Alimentarius and a quality management system that accommodates the requirements of EN ISO 9001 safeguard our exacting quality standards.

### Product quality

- Complies with the requirements of the European, American and Japanese pharmacopoeias
- Quality management system according to EN ISO 9001
- HACCP product safety system in accordance with the principles of the Codex Alimentarius
- Regular audits of the production facilities by experts and pharmaceutical companies
- Traceability system thanks to unique labelling

### Production process

esco HD-NaCl is obtained using a vacuum salt method from special and extremely pure brine. The sodium chloride obtained is of outstanding purity, contains no additives and is packaged separately from all other salt production equipment in dedicated areas. With this precise production process, we guarantee our customers maximum reliability and safety for products based on HD-NaCl.

### The company

esco – european salt company is the leading salt producer for international pharmaceutical companies in Europe and exports significant quantities overseas. With three state-of-the-art facilities in Europe, we are able to guarantee our customers reliable deliveries.